

## Product datasheet for **RR200332L4V**

### Cnksr2 (NM\_001113366) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Cnksr2 (NM_001113366) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Cnksr2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001113366
ORF Size:	3096 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200332).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001113366.1</a> , <a href="#">NP_001106837.1</a>
RefSeq Size:	3272 bp
RefSeq ORF:	3099 bp
Locus ID:	59322
UniProt ID:	<a href="#">Q9Z1T4</a>
Cytogenetics:	Xq21
Gene Summary:	membrane-associated guanylate kinase-interacting protein; may bind the postsynaptic density (PSD)-95 complex and synaptic scaffolding molecules [RGD, Feb 2006]



[View online »](#)