

Product datasheet for **RR200345L4V**

Kcnp2 (NM_001033961) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Kcnp2 (NM_001033961) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Kcnp2
Synonyms:	Kchip2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001033961
ORF Size:	660 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200345).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001033961.1 , NP_001029133.1
RefSeq Size:	2220 bp
RefSeq ORF:	663 bp
Locus ID:	56817
UniProt ID:	Q9JL59
Cytogenetics:	1q54
Gene Summary:	Ca(2+)-binding protein that acts as an auxiliary subunit for voltage-gated K(+) channels [RGD, Feb 2006]



[View online »](#)