

## Product datasheet for **RR200359**

### Cr1l (NM\_001005330) Rat Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Cr1l (NM_001005330) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Cr1l
Synonyms:	Cr1; Crry
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

**ORF Nucleotide Sequence:**

>RR200359 representing NM\_001005330  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGGAGGCTTCTCGCCTCTGGACCCCGTGGGCGCCTTGTAGCCTTCTGCCGCGGAGGAGTGCATCTGG  
 CGTCTGCTGCTGTTCTTGTGCGCATCTACTTTGGGCCAGTGCCAGCCCCACCACTGTTTCCTTATGC  
 CAAACCTATAAATCCAACGTATGAATCCACGTTTCCCGTTGGAACATCTTTGAAGTATGAATGTCGTCCA  
 GGATATATCAAGAGGCAGTTCTCTATCACCTGCGAAGTAACTCAGTGTGGACAAGTCTCAAGACGTGT  
 GTATACGTAACAATGTGAACTCCTTTAGATCCTCAGAATGGCATAGTACATGTAACACAGACATCCG  
 GTTTGGATCCTCTACTTATACCTGTAATGAAGGATACCGCCTCATTGGTTCTCCTCTGCTATGTGC  
 ATAATCTCTGATCAGAGCGTTGCCTGGGATGCCGAGGCACCTATTTGTAATCAATTCCTGTGAGATAC  
 CCCAAGCATTCCAATGGAGATTTCTCAGTCTAACAGAGAAGATTTTCATTATGGAATGGTAGTTAC  
 CTACCAGTGCAACTGATGCGAGAGGGAAGAAGCTCTTAACTGGTGGGTGAGCCCTCCATACACTGT  
 ACCAGCATCGATGGTCAAGTTGGAGTCTGGAGTGGCCCTCCTCCTCAGTGCATTGAACTCAACAAATGTA  
 CTCTCCCATGTTGAAAATGCAGTCATAGTGTCTAAAAACAAAAGCTGTTTTCTTAAGGGATATGGT  
 GGAGTTTATAGTGTGAGGATGGCTTTATGATGAAAGGAGACAGCAGTGTGTATTGTCGATCCCTAACAGA  
 TGGGAGCCTCAGTTACCAAGCTGCTTCAAGGTGAAATCCTGTGGTGTCTTCTGGGTGAGCTTCTAATG  
 GGCATGTATTCGTTCCACAAAATCTTCAACTCGGGGCCAAAGTGACCTTTGTTTGAATACAGGGTATCA  
 ATTAAAGGCAATTTCTTAGTCATTGTGTTCTTGATGGAGTGGAAAGCATTGGAATAGCAGCGTTCCT  
 GTGTGTGAACAAGTATGTAACCTCCCTCAGGATATGAGTGGATTCCAGAAAGGGTTGCAAAATGAAAA  
 AAGATTACTATGGAGATAATGTAGCCTTGGAAATGTGAGGATGGGTACTCTAGAAGGCAGTTCTCA  
 GAGCCAGTGCCAGTCCGATGCCAGCTGGATCCTCCTCTGCCCAAATGTGTCTCTCAAGTATATGTA  
 CTCCCTCAGGATATGAGTGGATTCCAGAAAGGGTTGCAAAATGAAAAAAGATTATTACTATGGAGATAATG  
 TAGCCTTGGAAATGTGAGGATGGGTATACTCTAGAAGGCAGTTCTCAGAGCCAGTCCAGTCCGATGCCAG  
 CTGGGATCCTCCTCTGCCCAAATGTGTCTCTCAAGTATGTAACCTCCCTCAGGATATGAGTGGATTC  
 CAGAAAGGGTTGCAAAATGAAAAAAGATTATTACTATGGAGATAATGTAGCCTTGGAAATGTGAGGATGGGT  
 TACTCTAGAAGGCAGTTCTCAGAGCCAGTCCAGTCCGATGCCAGCTGGGATCCTCCTCTGCCCAAATG  
 TGTCTCTCGTCAAACAGTGGTCTAATAGCTGGAATTTTCATTGGGATAATCGTCTTATTTTATTCATC  
 ATTTTTCTACTGGATGATTATGAAGTTTAAAAACGCAATCCACCAATGAAAAGTGTAAAGAAGTGG  
 GTATCTATTTAAATTCTAAAGAAGACAGCTGTGTTGAGCCTCAGTCTCTGCTCACAAGTCAGGAGAACAA  
 CAGTACCAGTAGCCAGCACGGAATCACTCACTCAAGAAGTC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:**


>RR200359 representing NM\_001005330  
 Red=Cloning site Green=Tags(s)

MEASSPLDPVGRVLVAFRCGGVHLAVLLFLSPSTLQPCPAPPLFPYAKPINPTDESTFPVGTSLKYECP  
 GYIKRQFSITCEVNSVWVSPQDVCIRKQCETPLDPQNGIVHVNTDIRFGSSITYTCNEYRILIGSSAMC  
 IISDQSVAWDAEAPICESIPCEIPPSIPNGDFFSPNREDFHYGMVVYQCNTDARGKKLFLNLVGEPSIHC  
 TSIDGQVGVWSGPPQCIENKCTPPHVENAVIVSKNLSFLSRDMVEFRCQDGFMMKGDSSVYCRSLNR  
 WEPQLPSCFKVKSCGAFGLGELPNGHVFPQNLQLGAKVTFVCNTGYQLKGNSSSHCVLDGVEIWNSSVP  
 VCEQVICKLPQDMSGFQKGLQMKKDYDDYDGNVALECEDGYTLEGSSQSQSDASWDPPLPKCVSQVICK  
 LPQDMSGFQKGLQMKKDYDDYDGNVALECEDGYTLEGSSQSQSDASWDPPLPKCVSQVICKLPQDMSGF  
 QKGLQMKKDYDDYDGNVALECEDGYTLEGSSQSQSDASWDPPLPKCVSRNSGLIAGIFIGIIVLILFI  
 IFSYWMIMKFKRNSTNEKCKEVIYLNKEDSCVQPQSLTSSQENNSTSSPARNSLTQEV


**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM\_001005330

ORF Size: 1863 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001005330.1</a> , <a href="#">NP_001005330.1</a>
<b>RefSeq Size:</b>	2280 bp
<b>RefSeq ORF:</b>	1866 bp
<b>Locus ID:</b>	54243
<b>UniProt ID:</b>	<a href="#">Q63135</a>
<b>Cytogenetics:</b>	13q27
<b>MW:</b>	68.6 kDa
<b>Gene Summary:</b>	This gene encodes a membrane-bound inhibitor of complement activation, which contains several tandemly arranged short consensus repeats (SCRs) of approximately 60 amino acids. The encoded protein may play a role in tumor immune surveillance, which is unrelated to complement inhibition. Three alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, Jul 2008]