

Product datasheet for **RR200390L4V**

Celf2 (NM_017197) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Celf2 (NM_017197) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Celf2
Synonyms:	Brunol3; Cugbp2; Etr3; Napor
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_017197
ORF Size:	1596 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200390).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_017197.3 , NP_058893.2
RefSeq Size:	2049 bp
RefSeq ORF:	1599 bp
Locus ID:	29428
UniProt ID:	Q792H5
Cytogenetics:	17q12.3
Gene Summary:	RNA binding protein that binds to two sets of AU-rich sequences in COX2; involved in mRNA stabilization and translational silencing [RGD, Feb 2006]



[View online »](#)