

Product datasheet for **RR200399L1V**

Pparg (NM_001145366) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pparg (NM_001145366) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Pparg
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_001145366
ORF Size:	1425 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200399).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001145366.1 , NP_001138838.1
RefSeq Size:	1830 bp
RefSeq ORF:	1428 bp
Locus ID:	25664
UniProt ID:	O88275
Cytogenetics:	4q42
Gene Summary:	ligand-activated transcription factor; mediates expression of genes involved in lipid metabolism [RGD, Feb 2006]



[View online »](#)