

Product datasheet for **RR200410L4V**

IIIrl1 (NM_001127689) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	IIIrl1
Synonyms:	Fit-1; FIT1; St2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001127689
ORF Size:	1698 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RR200410).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [NM_001127689.1](#), [NP_001121161.1](#)

RefSeq Size: 2065 bp

RefSeq ORF: 1701 bp

Locus ID: 25556

UniProt ID: [Q62611](#)

Cytogenetics: 9q22



Gene Summary: may be involved in modulating IL-1 signaling [RGD, Feb 2006]