

## Product datasheet for **RR200415L3V**

### Cav1 (NM\_133651) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Cav1 (NM_133651) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Cav1
Synonyms:	Cav
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_133651
ORF Size:	441 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200415).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_133651.2</a>
RefSeq Size:	2576 bp
RefSeq ORF:	444 bp
Locus ID:	25404
Cytogenetics:	4q21
Gene Summary:	structural protein of lipid raft domains in the membrane called caveolae; functions as a cholesterol-binding/shuttling protein [RGD, Feb 2006]



[View online »](#)