

Product datasheet for **RR200418L3V**

Cast (NM_001033716) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cast (NM_001033716) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Cast
Synonyms:	MGC108997
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001033716
ORF Size:	1956 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200418).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001033716.1 , NP_001028888.1
RefSeq Size:	4466 bp
RefSeq ORF:	1959 bp
Locus ID:	25403
UniProt ID:	P27321
Cytogenetics:	2q11
Gene Summary:	inhibitor of calpain; degradation is required for myoblast fusion; exists in multiple alternatively spliced isoforms [RGD, Feb 2006]



[View online »](#)