

Product datasheet for **RR200419L3V**

Cacna2d1 (NM_001110847) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cacna2d1 (NM_001110847) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Cacna2d1
Synonyms:	Cacna2; CCHLA2; DHSCCA
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001110847
ORF Size:	3237 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200419).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001110847.1 , NP_001104317.1
RefSeq Size:	3264 bp
RefSeq ORF:	3240 bp
Locus ID:	25399
Cytogenetics:	4q12
Gene Summary:	subunit of voltage-gated calcium channels; contributes to the regulation of N- and/or L-type calcium channels [RGD, Feb 2006]



[View online »](#)