

Product datasheet for **RR200425**

Adarb1 (NM_001111057) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Adarb1 (NM_001111057) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Adarb1
Synonyms:	Adar2; Red1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RR200425 representing NM_001111057
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCCCCCAAGGACAGCAGCACACCCGGGCTGGCGAGGGTATTCGGCTCTCCAACGGGGTGGTGGGA
 GCACCAGCAGGAAGCGGCCCTGGAGGAGGGCAGCAATGGCCACTCCAAGTACCGCTGAAGAAGCGAAG
 GAAAACGCCAGGGCCCGTTCTACCCAAGAACGCCCTGATGCAGCTGAACGAGATCAAACCGGCTTACAG
 TACATGCTGCTGTCCCAGACAGGACCCGTGCACGCACCTCTGTTTGTCTGTGTGGAGGTAAACGGGC
 AGGTCTTTGAAGGCTCCGGCCCTACAAAGAAGAAGGCAAAGCTGCATGCTGCTGAGAAGGCCCTGCGGTC
 TTTTGTCCAGTTTCCAACGCCTCTGAGGCCACCTGGCCATGGGAAGGACCTCTCCGTGAACACAGAC
 TTCACGTCGACAGGCGGACTTCCCGACACGCTCTTCAATGGCTTTGAGACTCCAGACAAGTCGGAGC
 CACCTTCTACGTAGGCTCCAATGGGGATGACTCCTTCAGCTCAAGCGGAGACGTTAGCCTGTCAGCCTC
 CCCAGTGCTGCCAGCCTTACCCAGCCTCTCTGCCATCCCACCACCATTCCCACCCCAAGTGGGAAG
 AACCCCGTGATGATCTTGAATGAGCTGCGCCAGGGCTGAAGTATGACTTCTCTCCGAGAGTGGGAGA
 GCCACGCCAAGAGCTTTGTCTGTCCGTGGTGGTATAGTGGCCAGTTCTTTGAGGGCTCAGGGAGAAACAA
 GAAGCTTGCCAAGGCCGGGCTGCACAGTCTGCCTTGGCTACTGTCTTCAATTTGCACTTGGACAAACG
 CCATCTCGCCAGCCTGTCTCAGTGAGGGTCTCCAGTTGCATTTGCCACAGGATTGGCAGATGCTGTCT
 CACGCCTGGTCTGGTAAGTTCAGTGACCTGACAGACAACCTTTCTCCCTCACGCACGAAGAAAAGT
 GCTCTCTGGAGTAGTGATGACCACAGGTACAGATGTCAAAGATGCCAAGGTGATAAGTGTTCGACAGGG
 ACGAAGTGCATCAACGGGAATACATGAGTGACCGTGGCTGGCTCTCAATGACTGCCACGCAGAGATAA
 TCTCCGAAGTCCCTGCTCAGTTTCTATACGCACAGCTCGAGCTTTACTTAAATAACAAGAAGACCA
 GAAAAAGTCCATATTTTCAGAAGTCAGAGCGGGTGGGTTCCGGCTGAAGGATACCGTGCAGTTCCACCTG
 TACATCAGCACCTCACCTGCGGAGACGCCAGAATATTCTCTCCCATGAGCCCGTGTAGAGGAACCGG
 CTGATAGACATCCGAATCGCAAAGCAAGGGGACAGCTGCGGACTAAAATAGAATCTGGCGAGGGGACAAT
 CCCTGTGCGCTCAAATGCCAGCATCCAGACCTGGGATGGGGTGTGCAGGGGGAACGGCTGCTCACCATG
 TCCTGCAGTGACAAGATAGCACGCTGGAACGTGGTGGCATCCAGGGGGCCCTGCTCAGCATTTTCGTGG
 AGCCCATCTACTTCTCCAGCATCATCTGGGCAGCCTGTACCACGGGGACCACCTCTCCAGGGCCATGTA
 CCAGCGGATCTCCAACATAGAGGACCTGCCACCCTCTACACCCTCAACAAGCCCTGCTCAGCGGTATC
 AGCAATGCAGAGGCACGGCAGCCAGGGAAGGCACCCAACCTCAGTGTCAACTGGACGGTGGGCGACACGG
 CCATTGAGGTCAATGCCACAACAGGGAAGGATGAGCTAGGCCGCCCTCCCGCTGTGTAAGCACGC
 GCTGTAAGTGTGCTGGATGCGGGTACACGAAAGGTGCCCCCACTGCTGCGCACCAAGATCACCAAG
 CCCACCACCTACCACGAGTCCAAGCTGGCAGCGAAGGAGTACCAGGCTGCCAAGGCACGCTGTCTACTG
 CCTTATCAAGGCGGGGCTGGGCGCTGGGTGGAGAAGCCACAGAGCAGGACCAGTTCTCTTCACTCC
 C

ACGCGTACGCGGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR200425 representing NM_001111057
Red=Cloning site Green=Tags(s)


MPPKDSSTPGPEGIPLSNGGGGSTSRKRPLEEGSNHSHKYRLKKRRKTPGPVLPKNALMQLNEIKPGLQ
 YMLLSQTGPVHAPL FVMSVEVNGQVFEVSGPTKKKAKLHAAEKALRSFVQFPNASEAHLAMGRTL SVNTD
 FTSDQADF PDTL FNGFETPKSEPPFVYVSGNGDDSFSSSGDVSL SASVVPASLTQPPLPIPPPPPPSGK
 NPVMILNELRPGLKYDFLSESGESHAKSFVMSVVVDGQFFEGSGRNKLLAKARAAQSALATVFNHLHDQT
 PSRQPVLSEGLQLHL PQVLADAVSRLVLGKFSDL TDNFSSPHARRKVL SGVMMTTGTDVKDAKVISVSTG
 TKCINGEYMSDRGLALNDCHAEIISRRSLLRFLYAQLELYLNNKEDQKKSIFQKSERGGFRLKDTVQFHL
 YISTSPCGDARIFSPHEPVLEEPADRHPNRKARGQLRTKIESGEGTIPVRSNASIQTDWGVLQGERLLTM
 SCSDKIARWNVVGIQGALLSIFVEPIYFSSIIILGSLYHGDHL SRAMYQRISNIEDLPPLYTLNKPLL SGI
 SNAEARQPGKAPNF SVNWTVGDTAIEVINATTGKDELGRPSRLCKHAL YCRWMRVHGKVPPLLRTKITK
 PTTYHESKLAKEYQAARLFTAFIKAGLGAWVEKPTEQDQFSFTP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



ACCN: NM_001111057

ORF Size: 2031 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001111057.1](#), [NP_001104527.1](#)

RefSeq Size: 6429 bp

RefSeq ORF: 2034 bp

Locus ID: 25367


UniProt ID: [P51400](#)

Cytogenetics: 20p12

MW: 74.1 kDa

Gene Summary: This gene encodes a double-stranded-RNA-specific adenosine deaminase that is involved in editing pre-mRNAs by site-specific conversion of adenosine (A) to inosine (I). Substrates for this enzyme include ionotropic glutamate receptors (GluR2-6) and serotonin receptor (5HT2C). Studies in rodents have shown that this protein can modify its own pre-mRNA by A->I editing to create a novel acceptor splice site, alternative splicing to which results in down regulation of its protein expression. Additional splicing events result in transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RR200425