

Product datasheet for **RR200443**

Bsg (NM_001109882) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Bsg (NM_001109882) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Bsg
Synonyms: 5A11; basigin; CE9; EMMPRIN; HT7; Ox47R
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR200443 representing NM_001109882
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGGATCGCC

ATGGCGGGCGGCTGCTGCTGGCGCTGGCCTTACGTTCTGAGTGGCCAAGGCGCTGCGGGCGGCTG
GCTTTCTCAAGGCGCCAATGTCCAGGAGCAGTGGCGGGGGCAGCGTGGTCTGCACTGTGAGGCTGT
GGGCAGCCCCATGCCGGAGATCCAGTGGTGGTTTGAAGGGAATGAACAAACGACAGCTGCTCCAGCTC
TGGGATGGCGCCCGGCTGGACCGTGTTCACATCCATGCCACCTACCGCCAGCATGCAGCCAGCACGCTGT
CTGTTGATGGGCTCGCCGAGAGGACACAGGCACTTATGAATGCCGGGCCAGCAGTGACCCTGACCGCAA
CCACCTGACCCGGCCACCCAGGGTCAAGTGGGTTCTGTCGCCAGGCGAGCGTGGTGGTCTTGAACCGGGC
ACCATCGTAACCTCTGTCCAGGAAGTTGACTCCAAGACACAGCTTACCTGCTTTTGAACAGCAGTGGCA
TTGACATCGTTGGCCACCCTGGATGAGAGGTGGCAAGTACTGCAGGAAGACACGCTGCCCGATCTACA
GATGAAGTACACGGTGGATGCAGATGACCGCTCTGGAGAATATCCTGCATCTTCCTTCTGAGCCTGTG
GGCAGAGGCAACATCAATGTGGAGGGGCCACCCAGGATCAAGTGGGAAAGAAATCGGAACACGCCAGTG
AGGGAGAGTTTGTGAAGCTGATCTGCAAGTCTGAGGCGTCCCACCCTCTGTGGATGAGTGGGTCTGGTT
TAAGACCTCTGACTGGGGACCAGACTATCTCCAACGGCACTGAGGCCAATAGCAAGTACGTATTATA
TCCACGCCTGAGCTGTGACAGCTGATCATCAGCGACCTGGACATGAACGTGGACCCTGGCACCTACGTGT
GCAATGCCACCAACTCCCAGGGCAGTGCTCGGGAGACCATCTCACTGCGTGTGCGCAGCCGCTGGCAGC
CCTCTGGCCCTTCTGGGCATTGTGGCCGAGGTCTGGTGTGGTACCATCATCTTCATCTACGAGAAG
AGGCGGAAGCCGGACCAGACCCTGGACGAGGATGATCCTGGCGCCGCCCACTGAAGGGCAGCGGTCTC
ACCTGAATGACAAGGACAAGAATGTGCCAGAGGAACGCCACC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR200443 representing NM_001109882
Red=Cloning site Green=Tags(s)

MAAALLLALAF TFLSGQGACAAAGFLKAPMSQEQWAGGSVVLHCEAVGSPMPEIQWWFEGNEPNDSCSQL
 WDGARLDRVHIHATYRQHAAS TLSVDGLAAEDTGT YECRASSDPDRNHL TRPPRVKQWVRAQASVVVLEPG
 TIVTSVQEVDSKTQLTCFLNSSGIDIVGHRWMRGGKVLQED TLPDLQMKYTV DADDRSGEYSCIFLPEPV
 GRGNINVEGPPRIKVGKKSEHASEGEFVKLICKSEASHPPVDEWVWFKTS DTGDQTI SNGTEANSKYVII
 STPELSELIISDLDMNVDPGT YVCNATNSQGSARETISLRVRSRLAALWPFLGIVAEVLVLTIIIF IYEK
 RRPKPDQLDEDDPGAAPLKGSGSHLNDKDKNVRQRNAT

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001109882

ORF Size:	1164 bp
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001109882.1 , NP_001103352.1
RefSeq Size:	1681 bp
RefSeq ORF:	1167 bp
Locus ID:	25246
UniProt ID:	P26453
Cytogenetics:	7q11
MW:	42.4 kDa
Gene Summary:	MRC OX-47 antigen that is upregulated on activated lymphocytes; member of the immunoglobulin superfamily [RGD, Feb 2006]