

Product datasheet for **RR200444L3V**

St6gal1 (NM_001113344) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	St6gal1 (NM_001113344) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	St6gal1
Synonyms:	Siat1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001113344
ORF Size:	1209 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200444).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001113344.1 , NP_001106815.1
RefSeq Size:	3879 bp
RefSeq ORF:	1212 bp
Locus ID:	25197
UniProt ID:	P13721
Cytogenetics:	11q23
Gene Summary:	may be involved in cell-type specific terminal glycosylation [RGD, Feb 2006]



[View online »](#)