

Product datasheet for **RR200447**

Fhl1 (NM_001033926) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Fhl1 (NM_001033926) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Fhl1
Synonyms: SLIM1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR200447 representing NM_001033926
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGCTTCTCAAAGACTCAGGTCCCTCCAGCTATAAGGTGGGCACCATGTCTGAGAAGTTCGACTGTC
ACTACTGCAGGGACCCCTTGCAGGGGAAGAAGTATGTGCAGAAGGATGGCCGTCCTGCTGCCTGAAGTG
CTTTGACAAGTTCTGCGCAACACCTGCGTGGAATGCCGCAAGCCATAAGCGCTGATGCTAAGGAGGTG
CACTATAAGAATCGTACTGGCAGCACCTGCTTCCGCTGTGCCAAGTGCCTTACCCCTTGGCCAGTG
AGACCTTTGTGTCCAAGGATGGCAAGATCCTGTGCAACAAGTGGCCACTCGGGAGGACTCCCCAGGTG
CAAAGGGTGCTTCAAGGCCATTGTGGCAGGAGACCAGAATGTGGAGTACAAGGGCACCATCTGGCATAAA
GACTGCTTACCTGCAGCAACTGCAAGCAAGTCATTGGGACCGGAAGCTTCTTCCGAAAGGGGAGGACT
TCTACTGTGTGACTTGCCATGAGACCAAGTTTCCCAAGCATTGCGTGAATGCAACAAGGCCATCACATC
TGGAGGAATCACTTACCAGGATCAGCCCTGGCATGCCGAGTGCTTTGTGTGTGTTACCTGCTCTAAGAAG
CTGGCGGGGCAGCGTTTACCGCTGTGGAGGACCAGTATTACTGCGTGGATTGCTACAAGAACTTTGTGG
CCAAGAAGTGTGCTGGATGCAAGAACCCCATCACTGGGTTTGGTAAAGGCTCCAGTGTGGTGGCCTATGA
AGGACAATCCTGGCAGCACTACTGCTTCCACTGCAAAAAATGCTCCGTAATCTGGCCAACAAGCGCTTT
GTCTTTCATAATGAGCAGGTGATTGTCCTGACTGTGCCAAAAAGCTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR200447 representing NM_001033926
 Red=Cloning site Green=Tags(s)

MASQRHSGPSSYKVGTMSEKFDCHYCRDPLQGKKYVQKDRHCCLKCFDKFCANTCVECCKPISADAKEV
 HYKNRYWHDTCFRCAKCLHPLASETFVSKDGIKCNKCATREDSRCKGCFKAIVAGDQNVYKGTIWHK
 DCFTCSNCKQVIGTGSFFPKGEDFYCVTCHETKFAKHCVKCNKAITSGGITYQDQPWHAECFVVCVTSKK
 LAGQRFTAVEDQYYCVDYKKNFVAKKACGCKNPITGFGKGSVVAYEGQSWHDYCFHCKKCSVNLANKRF
 VFHNEQVYCPDCAKKL


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001033926

ORF Size: 888 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001033926.2 , NP_001029098.1
RefSeq Size:	2333 bp
RefSeq ORF:	891 bp
Locus ID:	25177
Cytogenetics:	Xq37
MW:	33.6 kDa
Gene Summary:	contains a zinc binding domain; expressed in skeletal muscle [RGD, Feb 2006]