

## Product datasheet for **RR200451L4V**

### **Kng2 (Kng2l1) (NM\_012741) Rat Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Kng2 (Kng2l1) (NM_012741) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Kng2l1
Synonyms:	kininogen; KINKG; KINKH; Kng1; Kng2; Kngk
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_012741
ORF Size:	1917 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200451).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_012741.1</a> , <a href="#">NP_036873.1</a>
RefSeq Size:	2100 bp
RefSeq ORF:	1920 bp
Locus ID:	25087
UniProt ID:	<a href="#">P08934</a>
Cytogenetics:	11q23
Gene Summary:	precursor protein of kinin which is found in plasma; cysteine protease inhibitor and a major acute phase reactant [RGD, Feb 2006]



[View online »](#)