

Product datasheet for **RR200469**

Bcl2l1 (NM_031535) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Bcl2l1 (NM_031535) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Bcl2l1
Synonyms: bcl-X; Bcl-xl; Bcl2l; Bclx
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR200469 representing NM_031535
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGTCTCAGAGCAACCGGGAGCTGGTGGTTGACTTTCTCTCTACAAGCTCTCCCAGAAAGGATACAGCT
GGAGTCAGTTTAGCGATGTCGAAGAGAACAGGACTGAAGCCCCAGAAGAACTGAACCAGAAAGGGAGAC
CCCCAGTGCCATCAATGGCAACCCATCCTGGCACCTGGCGGATAGCCCCGCGGTGAATGGAGCCACTGGC
CACAGCAGCAGTTTGGATGCGCGGGAGGTAATCCCATGGCAGCAGTGAAGCAAGCGCTGAGAGAGGCTG
GCGATGAGTTTGAAGTACCGGATACCGGAGAGCATTCAAGTATCAACATCCCAGCTTCATATAACCCAGG
GACAGCATATCAGAGCTTTGAACAGGTAGTGAATGAACCTTTTCGGGATGGGGTAAACTGGGGTCGCATT
GTGGCCTTCTTCTCTTTGGCGGGCACTGTGCGTGAAAGCGTAGACAAGGAGATGCAGGTATTGGTGA
GTCGGATTGCAAGTTGGATGGCCACCTACCTGAATGACCACCTAGAGCCTTGGATCCAGGAGAACGGCGG
CTGGGTAAGAACCACGCCCTTGTGTGCCGCCCTTGTGTGTCTCTCTGTGGAGATCCCTAACTGC
CCTTTTGGTCTCTGCGCATGGTTGTTGAAGATATCGATTATTCAGGAGACATTCCTGGCTTCACTTTAA
TACCAGGGTTAACTTTGGGAATATTGATGACCCTGTTTTTAAGGAACCTGATTTTTTATTCTGGCTAC
CCTGTGGCCCCGAGTTTCATAGTTTTGTTCCAATTTCTCGGCAAGAAAAACAGCCTGTGTGTTACT
TGCTTAAAACC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR200469 representing NM_031535
 Red=Cloning site Green=Tags(s)

MSQSNRELVVDFLSYKLSQKGYWSQFSDVEENRTEAPEETEPERETPTSAINGNPSWHLADSPAVNGATG
 HSSSLDAREVIPMAAVKQALREAGDEFELRYRRAFSDL TSQLHITPGTAYQSFEQVYNELFRDGVNWGRI
 VAFFSFGGALCVESVDKEMQVLVSRIASWMATYLNHLEPWIQENGGWVRTTPLVCPPLVCLSSVEIPNC
 PFWSPGMVVEDIDYSGDIPGFTLIPGVNFGNIDDPVFKEPVFFILATLVAPQFHSPVPSRQRKTACVFT
 WLKT


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_031535

ORF Size: 852 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_031535.2 , NP_113723.2
RefSeq Size:	2203 bp
RefSeq ORF:	855 bp
Locus ID:	24888
UniProt ID:	P53563
Cytogenetics:	3q41
MW:	31.8 kDa
Gene Summary:	inhibits neuronal apoptosis; may provide neuroprotection against ischemic brain injury [RGD, Feb 2006]