

## Product datasheet for **RR200472L4V**

### Tpm1 (NM\_001034073) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Tpm1 (NM_001034073) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Tpm1
Synonyms:	Alpha-tm; Tma2; Tmsa
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001034073
ORF Size:	744 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200472).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001034073.1</a> , <a href="#">NP_001029245.1</a>
RefSeq Size:	1677 bp
RefSeq ORF:	747 bp
Locus ID:	24851
UniProt ID:	<a href="#">P04692</a>
Cytogenetics:	8q24



[View online »](#)

**Gene Summary:**

This gene is a member of the tropomyosin family and encodes a protein that heterodimerizes with a beta subunit, binding actin in both muscle and nonmuscle cells. In both smooth and striated muscle cells, it plays a role in calcium dependent-regulation of muscle contraction. In nonmuscle cells, it stabilizes cytoskeleton actin filaments. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008]