

Product datasheet for **RR200494L4V**

Pfkfb2 (NM_001033964) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pfkfb2 (NM_001033964) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Pfkfb2
Synonyms:	PFK-2; RH2K
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001033964
ORF Size:	1602 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200494).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001033964.2 , NP_001029136.1
RefSeq Size:	3755 bp
RefSeq ORF:	1605 bp
Locus ID:	24640
UniProt ID:	Q9JH5
Cytogenetics:	13q13
Gene Summary:	catalyzes the synthesis and degradation of fructose 2,6-bisphosphate [RGD, Feb 2006]



[View online »](#)