

Product datasheet for **RR200517L4V**

Eno1 (NM_001109908) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Eno1 (NM_001109908) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Eno1
Synonyms:	Nne
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001109908
ORF Size:	1302 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200517).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001109908.1 , NP_001103378.1
RefSeq Size:	1803 bp
RefSeq ORF:	1305 bp
Locus ID:	24333
UniProt ID:	P04764
Cytogenetics:	5q36



[View online »](#)

Gene Summary:

Multifunctional enzyme that, as well as its role in glycolysis, plays a part in various processes such as growth control, hypoxia tolerance and allergic responses. May also function in the intravascular and pericellular fibrinolytic system due to its ability to serve as a receptor and activator of plasminogen on the cell surface of several cell-types such as leukocytes and neurons. Stimulates immunoglobulin production.[UniProtKB/Swiss-Prot Function]