

Product datasheet for **RR200518L4V**

Crygc (NM_001081660) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Crygc (NM_001081660) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Crygc
Synonyms:	Cryg; Cryg3; Len
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001081660
ORF Size:	522 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200518).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001081660.1 , NP_001075129.1
RefSeq Size:	618 bp
RefSeq ORF:	525 bp
Locus ID:	24277
UniProt ID:	P02529
Cytogenetics:	9q32
Gene Summary:	putative protein of the lens of the eye; human homolog is lens protein likely involved in some forms of cataracts [RGD, Feb 2006]



[View online »](#)