

## Product datasheet for **RR200534L3V**

### **Itih3 (NM\_017351) Rat Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Itih3 (NM_017351) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Itih3
Synonyms:	PAIHC3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_017351
ORF Size:	2661 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200534).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_017351.1</a> , <a href="#">NP_059047.1</a>
RefSeq Size:	2783 bp
RefSeq ORF:	2664 bp
Locus ID:	50693
UniProt ID:	<a href="#">Q63416</a>
Cytogenetics:	16p16
Gene Summary:	heavy chain peptide component of pre-alpha inhibitor, a serum protein formed by covalent linkage of the chondroitin sulfate chain of bikunin to the heavy chain carboxyl-terminal residue [RGD, Feb 2006]



[View online »](#)