

## Product datasheet for **RR200564L3V**

### Raf1 (NM\_012639) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Raf1 (NM_012639) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Raf1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_012639
ORF Size:	1944 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200564).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_012639.2</a> , <a href="#">NP_036771.1</a>
RefSeq Size:	2970 bp
RefSeq ORF:	1947 bp
Locus ID:	24703
UniProt ID:	<a href="#">P11345</a>
Cytogenetics:	4q42
Gene Summary:	acts as a mitogenic protein kinase; mutant forms may play a role in transformation [RGD, Feb 2006]



[View online »](#)