

## Product datasheet for **RR200568L4V**

### Klf6 (NM\_031642) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Klf6 (NM_031642) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Klf6
Synonyms:	Aa1017; Copeb; CPBP; Zf9
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_031642
ORF Size:	951 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200568).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_031642.4</a> , <a href="#">NP_113830.2</a>
RefSeq Size:	4206 bp
RefSeq ORF:	954 bp
Locus ID:	58954
Cytogenetics:	17q12.1
Gene Summary:	activates transcription from the collagen alpha (I) promoter in vitro; may play a role in activation of hepatic stellate cells in response to injury [RGD, Feb 2006]



[View online »](#)