

Product datasheet for **RR200588**

Lmx1a (NM_001105967) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Lmx1a (NM_001105967) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Lmx1a
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR200588 representing NM_001105967
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTTGGACGGCCTGAAGATGGAAGAGAAGCTTTCAAAGTGGCATTGAGACCTCGGCCTCCTTCTCCTCTT
TGCTGGGCAGAGCGGTGAGCCCCAAGTCTGTCTGCGAGGGCTGTCAGCGGGTCATCTCGGATAGGTTTCT
GCTGAGACTCAACGACAGCTTCTGGCAGGCAATGCGTGCAGTGCCTCTTGCAAAGAGCCCCGGAG
ACCACCTGCTTCTACCGGACAAGAAGCTCTACTGCAAGTACCACTACGAGAAACTGTTTGTCTGCAAA
GTGGGGGTGCTTCAAGCCATTGCGCCCAATGAGTTTGTGTCATGCGTGCCAGAAGAGTGTGTACCACT
GAGCTGTTTCTGCTGCTGTCTGCGAGCGACAGCTGCAGAAGGGTGTGAGTTTGTCTGAAGGAGGGC
CAGCTGCTCTGCAAAGGAGACTATGAGAAGGAGCGAGAGTTGCTGAGCCTGGTGAAGCCAGCAGCTCAG
ACTCAGGCAAAAGTGTGATGAGGAGAGTCTTTGCAAGTCAAGCCACGGGGCAGGAAAGGAGCATCAGA
GGACGGCAAGGACCATAAGCGACCCAAACGTCCCAGAACCATCTTGACCACTCAGCAGAGGAGAGCATT
AAGGCCTCATTTGAAGTATCCTCAAGCCCTGCAGAAAGGTGAGGGAGACTCTGGCTGCGGAGACAGGGC
TGAGTGTCCGTGTGGTTCAGGTGTGGTCCAGAACCAGCGAGCCAAGATGAAGAAGCTGGCCCGGCGACA
GCAGCAACAGCAACAGGATCAACAGAACCAGAGGCTGAGTTCTGCTCAGACAAATGGTGTGGCAAT
GCTGGCATGGAAGGCATCATGAACCCCTATACAACGTTGCCACCCACAGCAACTGCTGGCCATTGAAC
AGAGCGTCTACAACCTGATCCCTTCCGACAGGGTCTACCCACCCAGATGCCTGGAGATCACATGCA
CCCCTACGGTGTGAACCTCTTTTCCATGACCTGGATAGTGACGACACATCTCTCAGTAACCTGGGGCAG
TGCTTCTGGCAACCTCAGAAGCTGGACCCCTACAGTCCAGAGTGGGAAACCCATCGACCATCTGTATT
CCATGCAGAATTCTATTACCTCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR200588 representing NM_001105967
 Red=Cloning site Green=Tags(s)


MLDGLKMEENFQSAIETSASFSSLLGRAVSPKSVCEGCQRVISDRFLLRLNDSFWHEQCVQCASCKEPLE
 TTCFYRDKKLYCKYHYEKLFAVKCGGCFEAIAPNEFVMRAQKSVYHLSCFCCCVCERLQKGFVFLKEG
 QLLCKGDYEKERELLSLVSPAASDSGKSDDEESLCKSAHGAGKGASEDGKDHKRPKRPRTILTTQRRFA
 KASFEVSSKPCRKVRETLAAETGLSVRVVQVWFQNRQAKMKKLARRQQQQDQDQNTQRLSSAQTNGSGN
 AGMEGIMNPYTTLPQPQLLAIIEQSVYNSDPFRQGLTPPQMPGDHMPYGAELFHDLDSDDTSLSNLGD
 CFLATSEAGPLQSRVGNPIDHLYSMQNSYFTS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN:

NM_001105967

| | |
|-------------------------------|---|
| ORF Size: | 1146 bp |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | NM_001105967.2 , NP_001099437.2 |
| RefSeq Size: | 3252 bp |
| RefSeq ORF: | 1149 bp |
| Locus ID: | 289201 |
| Cytogenetics: | 13q24 |
| MW: | 42.8 kDa |