

Product datasheet for **RR200596**

Mrpl19 (NM_001029898) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Mrpl19 (NM_001029898) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Mrpl19
Synonyms: MRP-L15; Rpml15
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR200596 representing NM_001029898
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCGGCCTCCATGGCAGAAAGTTGTAGGGTGTCCATGTACCTGGCCTGGAGCTACCGAATGGCCAGGC
CACGGCTTGCTGCTCTCGCTCTGATGCCTGCCGTGCCACACGGGCCCAAGTCGGTTTCAGAGCACTGG
GCCATCCGAGCCCGGTGGTTCAAGCCGCCGCGAAACCTGTTCATTGTGGACCGCGCCGAGCCCCAGAG
GAAGAGAGAAGGTTCTGAGTCCTGAATTCATCCCTCCAAGGGGAGAAACAAACCTCTGAAATTTGAGA
TAGAAAGAAGAGATATGCTGGACAGGAGAAAGGTGCTTCCATTCCAGAGTTCTATGTCCGGAAGCATACT
GCGCGTCAACCACAGCTGACCCCTATGCCAGTGGGAAAACCCAGCCAGTTTCTGGGGATTTGCATCAAGCGA
TCAGGAAGTGGGCTTGGGGCTACTTTCACCCTCAGGAACACGATCGAAGGACAAGGTGTGGAGATTTGCT
TTGAACTGTATAACCTCGAATCCATGAGATCCAAGTGGTGAAGTTGGAGAAGCGGCTGGATGACAGCCT
GCTGTACTTACGGGATGCCCTCCCCGAGTACAGCACCTTCGATGTGAATATGAAGCCAGTTCCACTTGAG
TCCTGCCAGGAAGTGCCTGTGAACAAGCTGAAAGTTAAAATGAAGCCAAAGCCGTGGTCTAAACGCTGGG
AACGTCCGAACCTCAACATTAAGGGATCAGATTTGATCTTAACTGAAGACAAAATGAAAGAAGC
CCAGAAGTGGAGCAAGCCATGGCTTGATTTGATATGATGAGAGAATATGACACTTCAAAAATTTGAAGCT
GCATTGTGGGAAGAAATAGAAGCATCAAAGAAGTCT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR200596 representing NM_001029898
 Red=Cloning site Green=Tags(s)

MAASMAESCRVSMYLAWSYRMAPRPLAALASDACRAHTGPSRFQSTGPSEPGGFKPPPVPVIVDRRRAPE
 EERRFLSPEFIPPRGRTNPLKQIERRDMLDRKVLPIPEFYVGSILRVTTADPYASGKTSQFLGICIKR
 SSGSLGATFTLRNTIEGQVEICFELYNPRIHEIQVVKLEKRLDSSLYLRLDALPEYSTFDVNMKPVPLE
 SCQEVPNKLVKMKPKPWSKRWERPNFNIGIRFDLALTEEQMKEAQKWSKPWLEFDMMREYDTSKIEA
 ALWEEIEASKKS


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001029898

ORF Size: 876 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001029898.1 , NP_001025069.1
RefSeq Size:	1169 bp
RefSeq ORF:	879 bp
Locus ID:	297372
Cytogenetics:	4q34
MW:	33.6 kDa
Gene Summary:	may act as a structural component of mitochondrial ribosomes [RGD, Feb 2006]