

Product datasheet for **RR200609L3V**

Acot7 (NM_013214) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Acot7 (NM_013214) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Acot7
Synonyms:	Bach
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_013214
ORF Size:	1116 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200609).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_013214.4 , NP_037346.2
RefSeq Size:	1453 bp
RefSeq ORF:	1119 bp
Locus ID:	26759
UniProt ID:	Q64559
Cytogenetics:	5q36
Gene Summary:	cleaves acyl-CoA thioesters to free fatty acids and coenzyme A [RGD, Feb 2006]



[View online »](#)