

## Product datasheet for **RR200618**

### Tnfsf4 (NM\_053552) Rat Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Tnfsf4 (NM\_053552) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Tnfsf4  
**Synonyms:** Ox40l; Tnlg2b  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR200618 representing NM\_053552  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGC**C

ATGGAAGGGGAAGGGTTCAACCCCGGATGAGAATCTGGAAAATGGATCAAGGCCAAGGTTCAAGTGA  
 AGAAGGTGCTAAGCTGGTGGTCTCTGGGATTAAGGCAGCAGGGCTGCTCCTGTGCGTCGTCTATGTCTG  
 CCTGCAATTCTTTCATCTCCGGCGAAGGACTCTCCAATCCAAAGACTCAGAGCACCAGTTACCGGATGT  
 GAGGGGGGAAGACTATTCAATGGCACATCCAAGAATGAGTATGAAACCATGGAGGTGCAGAACAAATTCGG  
 TCATCATCAACTGCGATGGGCTTTATCTCATCCACCTGAAGGGCTCCTTTTTCCAGGAGGTCAAGATCAA  
 CCTTCAATTTCCGGAAGGATCGCAGTCCAATCTTTGTGCCAATGCTGAACAATGGTCAAAGGGTTGTTTTT  
 ACTGTGGTGACCTCTTTGGCTTTTAAAGACGAAGTTACTTGACTGTAATGCTTCTGATACTCTCTGTG  
 AACACCTCCAGATAAATGATGGGGAAGTATTATGTCCAGCTAACTCCTAATGGATACTGTGCTCCTGA  
 AAGACCTTATAGCAGCACTGTGAACCAAGTACCACTG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RR200618 representing NM\_053552  
 Red=Cloning site Green=Tags(s)

MEGEGVQPPDENLENGSRPRFKWKVLRLLVSGIKAAGLLL CVVYVCLQFSSSPAKDSPIQRLRAPVTGC  
 EGGRLF IGTSKNEYETMEVQNNSVI INCDGL YL IHLKGSFFQEVK INLHFRKDRSPI FVPM LNNGQRVVF  
 TVVTSLAFKDEVYLT VNASDTL CEHLQ INDGEL IIVQLTPNGYCAPERPYSSTVNQVPL


**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI



[View online »](#)

**Cloning Scheme:**

**Plasmid Map:**


ACCN: NM\_053552

ORF Size: 597 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_053552.1](#), [NP\\_446004.1](#)

**RefSeq Size:** 1804 bp

**RefSeq ORF:** 600 bp

**Locus ID:** 89814

**UniProt ID:** [Q9Z2P3](#)

**Cytogenetics:** 13q22

**MW:** 22.2 kDa

**Gene Summary:** ligand for the tumor necrosis factor receptor OX40 (CD134); exhibits costimulatory activity for T-cell proliferation [RGD, Feb 2006]