

## Product datasheet for **RR200658L3V**

### Phf5a (NM\_138888) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Phf5a (NM_138888) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Phf5a
Synonyms:	Ac2-246
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_138888
ORF Size:	330 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200658).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_138888.1</a> , <a href="#">NP_620243.1</a>
RefSeq Size:	872 bp
RefSeq ORF:	333 bp
Locus ID:	192246
UniProt ID:	<a href="#">P83871</a>
Cytogenetics:	7q34
Gene Summary:	zinc-finger transcription factor; promotes the up-regulation of the connexin43 gene by estrogen [RGD, Feb 2006]



[View online »](#)