

## Product datasheet for **RR200789**

### S1pr2 (NM\_017192) Rat Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** S1pr2 (NM\_017192) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** S1pr2  
**Synonyms:** Edg5; Gpcr13; GPCR18; H218; snGPCR18  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR200789 representing NM\_017192  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGGGCGGTTTACTCAGAGTACCTCAATCCTGAGAAGGTTCAGGAACACTACAATTACACCAAGGAGA  
CGCTGGACATGCAGGAGACGCCCTCCCGAAGGTGGCCTCCGCCCTCATCATATTTTATGCTGTGCCAT  
CGTGGTGGAGAACCTTCTGGTGTAATCGCAGTGGCCAGGAACAGCAAGTCCACTCAGCCATGTACCTG  
TTCTCGGCAACCTGGCAGCCTCCGACCTGCTGGCAGGCGTGGCCTTCGTGGCCAACACCTTGCTCTCCG  
GACCTGTACCCTGTCCCTAACTCCCTTGACAGTGGTTTGGCCGAGAGGGTTCAGCCTTCATCAGCTCTC  
TGCTCGGTCTTCAGCCTCCTGGCCATTGCCATCGAGAGACAAGTGGCCATCGCCAAGGTCAAGCTCTAC  
GGCAGTGACAAAAGCTGTGCAATGTTGATGCTCATTGGGGCCTTGGCTGATATCGCTGATTCTGGGTG  
GCTTGCCCATCCTGGGCTGGAATTGTCTGGACCATCTGGAGGCTTGCTCCACTGTGCTGCCCTCTATGC  
TAAGCACTATGTGCTCTGCGTGGTCACCATCTTCTGTCTACTGGCTATCGTGGCCTGTACGTC  
CGAATCTACTTCGTAGTCCGCTCAAGCCATGCGGACGTTGCTGGTCTCAGACGCTGGCCTGCTCAAGA  
CAGTCACCATCGTACTGGGTGTTTTATCATCATCTGCTGGCTGCCGGCTTTAGCATCCTTCTTACTGACTC  
TACCTGTCCCCTCGGGCCTGTCTGTCTTACAAAGCCATTATTTCTTGGCCTCGCCACCCTCAAC  
TCTCTGCTCAACCCTGTCTATACATATGCGTAGCCGGACCTTCGGAGGGAGGTACTGAGGCCCTGC  
TGTGCTGGCGGCAGGGGAAGGGAGCAACAGGGCGCAGAGGTGGGAACCCTGGTACCAGCTCTGCCCT  
CCGCAGCTCCAGCTCCCTGGAGAGAGGCTTGATATGCCTACATCGCCAACATTTCTGGAGGGCAACACA  
GTGGTC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR200789 representing NM\_017192  
 Red=Cloning site Green=Tags(s)

MGGLYSEYLNPEKVQEHYNYTKETLDMQETPSRKVASAFIIILCCAIVVENLLVLIAVARNKSFHSAMYL  
 FLGNLAASDLLAGVAFVANTLLSGPVTLSTPLQWFAREGSAFITLSASVFSLLAIAIERQVAIAKVKLY  
 GSDKSCRMLMIGASWLSLILGGLPILGWNCLDHLEACSTVLPYAKHYVLCVVTIFSVILLAIVALYV  
 RIYFVVRSSHADVAGPQTLALLKTVTIVLGVFIIICWLPAFSILLLDSTCPVRACPVLKHAHYFFAFATLN  
 SLLNPVIYTWRSRDLRREVLRLLLCWRQKGATGRRGGNPGHRLPLRSSSSLERGLHMPTSPTFLEGNT  
 VV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN:

NM\_017192

<b>ORF Size:</b>	1056 bp
<b>OTI Disclaimer:</b>	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <a href="mailto:custsupport@origene.com">custsupport@origene.com</a> or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a></p>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_017192.1</a> , <a href="#">NP_058888.1</a>
<b>RefSeq Size:</b>	2754 bp
<b>RefSeq ORF:</b>	1059 bp
<b>Locus ID:</b>	29415
<b>UniProt ID:</b>	<a href="#">P47752</a>
<b>Cytogenetics:</b>	8q13
<b>MW:</b>	38.7 kDa
<b>Gene Summary:</b>	G protein-coupled receptor for sphingosine 1-phosphate [RGD, Feb 2006]