

## Product datasheet for RR200817L4

### Tmed3 (NM\_001004249) Rat Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Tmed3 (NM_001004249) Rat Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Tmed3
Synonyms:	MGC94283
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200817).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.



[View online »](#)



Locus ID: 300888

UniProt ID: [Q6AY25](#)

Cytogenetics: 8q31

Gene Summary: Potential role in vesicular protein trafficking, mainly in the early secretory pathway. Contributes to the coupled localization of TMED2 and TMED10 in the cis-Golgi network (By similarity).[UniProtKB/Swiss-Prot Function]