

Product datasheet for **RR200818L3V**

Gstm1 (NM_017014) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Gstm1 (NM_017014) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Gstm1
Synonyms:	GSTA3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_017014
ORF Size:	654 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200818).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_017014.1 , NP_058710.1
RefSeq Size:	1055 bp
RefSeq ORF:	657 bp
Locus ID:	24423
UniProt ID:	P04905
Cytogenetics:	2q34
Gene Summary:	multifunctional enzyme thought to have scavenger function; able to neutralize harmful metabolites from foreign substances [RGD, Feb 2006]



[View online »](#)