

Product datasheet for **RR200899L3V**

Sart1 (NM_031596) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Sart1 (NM_031596) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Sart1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_031596
ORF Size:	2418 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200899).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_031596.1 , NP_113784.1
RefSeq Size:	2532 bp
RefSeq ORF:	2421 bp
Locus ID:	29678
UniProt ID:	Q5XIW8
Cytogenetics:	1q43
Gene Summary:	gene expressed in the nucleus and upregulated in proliferating cells, also found in cytoplasm of epithelial cancer cells [RGD, Feb 2006]



[View online »](#)