

Product datasheet for **RR200909L4V**

Tra2b (NM_057119) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tra2b (NM_057119) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Tra2b
Synonyms:	Sfrs10; Srsf10
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_057119
ORF Size:	864 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200909).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_057119.1 , NP_476460.1
RefSeq Size:	1978 bp
RefSeq ORF:	867 bp
Locus ID:	117259
UniProt ID:	P62997
Cytogenetics:	11q23
Gene Summary:	may act as an RNA splicing factor; may play a role in response to reactive oxygen species [RGD, Feb 2006]



[View online »](#)