

Product datasheet for **RR200911**

Nuak1 (NM_001106774) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Nuak1 (NM_001106774) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Nuak1
Synonyms:	RGD1309956
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide
Sequence:

>RR200911 representing NM_001106774
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGAAGGGCGGCAGTGTCCGCGGCAGGGGATGGCCCTGCTGTGGAGACGGGTTTGCCCGGCTCTCCGC
GGGAGGTAGTGGCCGGGCGACTGCTGCTCCAGTGAACCCAGGAAGCCGCACGGGTGAAGCGACATCA
CCACAAACACAACCTTGAAGCATCGCTATGAGCTGCAAGAGACCCTGGGCAAAGGTACCTACGGCAAAGTC
AAGAGAGCCACGGAGAGGTTTTTCAGGCCGAGTGGTTGCTATAAAATCCATCCGTAAGGACAAAATCAAGG
ATGAGCTAGACATGGTTCACATCAGACGAGAGATTGAGATTATGTCATCCCTCAACCATCCTCATATCAT
CAGTATTTATGAAGTGGTTGAGAACAAAGATAAGATTGTGATCATCATGGAATATGCCAGCAAAGGAGAG
CTGTATGACTACATCAGTGAGCGGAGACGCCTCAGCGAGAGGGAGACCAGGCACCTTCCGGCAGATCG
TGCTGCCGTGCATTACTGTCACAAGAATGGCGTGGTCCATCGGGACTTGAAGCTGGAGAACATACTGCT
TGATGACAACGTAAATATTAAGATCGCTGACTTTGGGCTTCCAACCTGTACCAGAAGGACAAGTCTTG
CAAACATTTTGTGGGAGCCCACTCTATGCATCACCCGAGATTGTCAATGGCAGGCCCTTACCGGGGGCCAG
AGGTGGACAGCTGGGCCCTGGGCGTCTTGTCTTACACTCTCATTTATGGAACAATGCCCTTCGATGGCTT
TGATCACAAAAACCTCATCCGACAGATCAGCAGTGGAGAGTACCGGGAGCCACACAACCCTCAGATGCT
CGGGGACTCATCCGGTGGATGCTGATGGTGAATCCTGACCGCCGGGCCACCATCGAGGACATCGCCAACC
ACTGGTGGGTGAAGTGGGCTACAAGAGCAGCGTCTGTGACTGCGATGCCCTTCCCGACTCCGAGTCTCC
GCTCTTGGCTCGGATCATTGATTGGCATCACCGCTCCACCGGGCTGCAGGCTGAGGCTGACGCCAAGATG
AAGCCAAGCCCCGGGCCCTCGGAGTTGTGTTGGAGCGCAACGGTCACTGAAGAAGTCCAAGAAAGAGA
ACGACTTTCCCAAGTCTGGCCAGGACTCATTGCCTGAAAGCCCATCCAAGCTGAGTCCAAGAGGCCCAA
GGGCATCCTGAAGAAGCGGAGCAACAGTGAACATCGCTCCCACAGCACTGGCTTCATCGAAGGCATCGTG
GGCCCCGCTTACCATCCCTTTCAAATGGAGCAAGACTTGTGCCGGACTGCCATCCCCCTCCCGGCT
CCCCGAGGCGCAGGGAAGCTCAGGGAAGCTCAGCCTCAAACAGTCTGCCACCATGCCAAAGAAAGGCAT
CTTGAAAAGACCCAGCAGAGAGAATCTGGTTACTACTCGTCCCCAGAGCGCAGCGAGTCTCAGAGCTG
CTGGACAGTAACAATGTGATCATCAGCAGTGGCTCTCTCCCCCTCCCCAGACCCAGCCAGGGGAACGT
CCCACAGCCTGCTCCTGCCGGAGGAAGGGCATCTTGAAACACAGCAGCAGGTAAGTCCGGATGGTGCCACAGA
CCCTGCCCTCACCGCCCCGAAATACCCACACTGGAATCCTTGTGCTCCCGGTGTCCTCCGATGGC
ATCTCTCGGAGCTACAGTCGCCCCGTCAGTATCATCAGTGACGACAGCGTCTGTCCAGCGACTCTTTCG
ACTTGTTAGATCTCCAGGAGAATCGGCCGCCCGCCAGCGCATCCGCAGCTGCGTCTCTGCGGAAAACCT
CCTCCAGTTCAGGACTTTGAGACTCCCCACAACCGGCCCGGCCCAAGTACCTGAAGCGGTACCGGAAC
CGGCTGGCAGACAGTAGCTTCTCCCTTCTGCTGATATGGACGACGTGACTCAGGTGATAAGAAGGCTC
TGGAGATCTGCAGCAAGCTCAAC


ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA


Protein Sequence: >RR200911 representing NM_001106774
 Red=Cloning site Green=Tags(s)

MEGAAVSAAGDPAVETGLPGSPREVVAGATAAPVEPRKPHGVKRHHHKKHNLKHRYELQETLGKGTYGKV
 KRATERFSGRVVAIKSIRKDKIKDELDMVHIRREIEIMSSLNPHIISIYEVFENKDKIVIIMEYASKGE
 LYDYISERRRLSERETRHFFRQIVSAVHYCHKNGVVHRDLKLENILLDDNCNIKIADFGLSNLYQKDKFL
 QTFCGSPLYASPEIVNGRPYRGPEVDSWALGVLLYTLIYGTMPFDGFDHKNLIRQISSGEYREPTQPSDA
 RGLIRWMLMVNPDRRATIEDIANHWWVNWGYKSSVDCDALPDSEPLLARIIDWHHRSTGLQAEADAKM
 KAKPGASEVVLERQRSLKSKKENDFPQSGQDSLPEPSKLSKRPKGILKKRSNSEHRSHSTGFIEGIV
 GPALPSPFKMEQDLCRTAIPPLGSPAQSSGKLSLKQSATMPKKGILKKTQQRESGYSSPERSESSEL
 LDSNNVISSGLSSPPDPARGTSHSLSCRRKGILKHSSRYSDGATDPALTGPEIPTLESLSPPGVPDGG
 ISRSYSRPSSIISSDVLSSDSFDLLDLQENRPARQRIRSCVSAENFLQLQDFETPHNRPQYLKRYRN
 RLADSSFLLADMDDVTQVYKKALEICSKLN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:
Cloning Scheme:



Plasmid Map:


ACCN: NM_001106774

ORF Size: 1983 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001106774.1](#), [NP_001100244.1](#)

RefSeq Size: 2864 bp

RefSeq ORF: 1986 bp

Locus ID: 299694

Cytogenetics: 7q13

MW: 74 kDa