

## Product datasheet for **RR200919L3V**

### Agfg1 (NM\_001135596) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Agfg1 (NM_001135596) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Agfg1
Synonyms:	MGC116375
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001135596
ORF Size:	1683 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200919).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001135596.1</a> , <a href="#">NP_001129068.1</a>
RefSeq Size:	3267 bp
RefSeq ORF:	1686 bp
Locus ID:	363266
UniProt ID:	<a href="#">Q4KLH5</a>
Cytogenetics:	9q34
Gene Summary:	Required for vesicle docking or fusion during acrosome biogenesis. May play a role in RNA trafficking or localization (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)