

## Product datasheet for RR200955L4

### Smpdl3b (NM\_001025737) Rat Tagged Lenti ORF Clone

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | Expression Plasmids  |
| Product Name:             | Smpdl3b (NM_001025737) Rat Tagged Lenti ORF Clone              |
| Tag:                      | mGFP   |
| Symbol:                   | Smpdl3b  |
| Synonyms:                 | MGC116117  |
| Mammalian Cell Selection: | Puromycin  |
| Vector:                   | pLenti-C-mGFP-P2A-Puro (PS100093)                              |
| E. coli Selection:        | Chloramphenicol (34 ug/mL)                                     |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(RR200955). |
| Restriction Sites:        | SgfI-MluI  |
| Cloning Scheme:           |  |

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.



[View online »](#)



Locus ID: 362619

Cytogenetics: 5q36