

## Product datasheet for **RR200957**

### **Srpra (NM\_001034150) Rat Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Srpra (NM_001034150) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Srpra
Synonyms:	Srpr
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide  
Sequence:

>RR200957 representing NM\_001034150  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCCGGATCGCC

ATGCTTGACTTCTTACCATCTTCTCGAAGGGCGGGCTTGTGCTCTGGTGTTCAGGGCGTGAGCGACT  
CGTGCACCGGGCCCGTGAACGCGTTGATTTCGTTCTGTACTGCTGCAGGAACGGGGAGGTAAACAACCTTT  
CACTCATGAGGCTCTCACACTCAAGTATAAACTGGACAACCGATTGAACTGGTGTGGTTGGCTTT  
CAGAAGATCCTCACACTGACCTATGTAGACAAATTGATAGATGATGTGCACCGGCTGTTCCGGGATAAGT  
ACCGCACAGAGATCCAACAACAGAGTGCCTTAAGTTTATTGAATGGCACTTTTGTATTTCCAAAATGACTT  
CCTGCGGCTCCTTCGTGAAGCAGAGGAGAGCAGTAAGATCCGTGCTCCCACTACCATGAAGAAATTTGAA  
GATTCTGAAAAAGCTAAGAACTGTGAGATCCATGATTGAAACACGGGGGAAAAGACCAAGGAAAAAG  
CAAAAAACAACAAAAAAGGGTCCAAAAAGGAAGTTCTGATGGCCCTTTGGCTACCAGCAAAACAGC  
CCCTGCAGAAAAGTCAGGTCTCCCACTGGGACCTGAGAATGGGGAACCTTCCAAAGAGGAGCTGATACGT  
AGGAAGAGAGAGGAGTTCATTCAGAAGCATGGGAAGGGCTGGACAAATCCAGTAAATCTACAAAGTCAG  
ATATTCAAAAGGAGAAAAGGCAAAAAAGCACCCCGGGTGTGGAACTAGGTGGCTGTGCTAACAAGGAAGT  
CTTGGATTACAGTACTCCACCACCAATGGAATCCAGAGGCAGCCCTTTCTGAAGATATCAACTTGATT  
CGAGGAACTGGGCTGGGGACAGCTTCAAGATCTGGATTGCAGCAGTTCAGATGATGAAGGGGCCACTC  
AAAACACCAAACTAGTCTACCAAGGGAACCTGCGGTGGCATGTTTGGGATGCTGAAGGGTCTGGTGGG  
TTCCAAGAGCTGAGTCGAGAAGACATGGAATCTGTGCTAGACAAGATGCGTGATCATCTCATTGCTAAG  
AATGTGGCTGCTGACATTGCTGTCCAGCTCTGTGAATCTGTTGCCAACAAGTTAGAAGGCAAAGTATGG  
GGACATTTAGCACTGTGACTTCAACAGTGAAGCAAGCTCTACAAGAGTCTCTGGTGCAGATCCTGCAGCC  
ACAGCGTCGAGTAGACATGCTCCGGGATATCATGGATGCCAAGCTCGCCAGCGCCCTTACGTAGTTACC  
TTCTGTGGTGTTAATGGAGTGGGGAAGTCTACTAATCTTGCCAAGATCTCCTTCTGGCTTTTAGAGAATG  
GCTTTAGCGTCTCATTGCTGCCTGCGATACATTTCTGCTGGGGCTGTTGAGCAGCTCCGCACACACAC  
CCGGCGCTGACTGCCCTTATCCCCGGAGAAGCATGGTGGCCGAACGATGGTGCAGTTGTTTAAAAA  
GGCTATGGCAAGGACGCTGCTGGCATTGCCATGGAAGCCATTGCCTTTCACGTAACCAAGGCTTTGACG  
TGCTGCTGGTGGACACTGCTGGCCGATGCAAGACAATGCTCCTCTGATGACTGCCCTGGCCAACTCAT  
TACTGTCAATACACCTGACTTGGTGTGTTTGTGGGGAGGCCTTAGTAGGCAATGAAGCTGTGGATCAA  
CTGGTCAAGTTCAACAGAGCCTTAGCTGACCATTCTATGGCTCAGACACCTCGGCTCATTGATGGCATTG  
TTCTAACCAAAATTTGACACCATTGATGACAAGTGGGAGCTGCTATTTCTATGACATACATCACAAGCAA  
ACCCATCGTCTTTGTGGTACTGGCCAGACCTACTGTGACCTACGCAGTCTCAATGCCAAGGCTGTGGT  
GCTGCCCTCATGAAGCT

AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC  
TGGATTACAAGGATGACGACGATAAGGTTTAA


Protein Sequence: >RR200957 representing NM\_001034150  
 Red=Cloning site Green=Tags(s)

MLDFFTIFSKGGLVWCFQGVSDSCTGPVNALIRSVLLQERGGNNSFTHEALTLKYKLDNQFELVFVVG  
 QKILTLTYVDKLIIDVHRLFRDKYRTEIQQSALSLNGTFDFQNDFLRLREAEESKIRAPTTMKKFE  
 DSEKAKKPVRSMIETRGEKTEKAKNNKKKGSKKEGSDGPLATSKTAPAEKSGLPVGPENGELSKEELIR  
 RKREEFIQKHGKGLDKSSKSTKSDIPKEKGKKAPRVWELGGCANKEVLDYSTPTTNGTPEAALSEDINLI  
 RGTGPGGQLQDLDCSSSDDEGATQNTKPSATKGTGGMFGMLKGLVGSLSREDMESVLDKMRDHLIAK  
 NVAADIAVQLCESVANKLEGKVMGTFSTVTSTVKQALQESLVQILQPQRRVMDLRDMDAQRQRYPVVT  
 FCGVNGVGKSTNLAKISFWLLENGFSVLIACDTFRAGAVEQLRTHTRRLTALHPPEKHGGRTMVQLFEK  
 GYGKDAAGIAMEAIAFARNQGFVVLDVDTAGRMQDNAPLMTALAKLITVNTPLVLFVGEALVGNEAVDQ  
 LVKFNRLADHSMAQTPRLIDGIVLTKFDTIDDKVGAAISMTYITSKPIVFGTGQTYCDLRSNAKAVV  
 AALMKA


SGPTRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-RsrII

Cloning Scheme:



Plasmid Map:



ACCN: NM\_001034150

ORF Size: 1908 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001034150.1](#), [NP\\_001029322.1](#)

RefSeq Size: 3001 bp  
RefSeq ORF: 1911 bp  
Locus ID: 315548  
Cytogenetics: 8q21  
MW: 69.7 kDa