

## Product datasheet for **RR200995L4V**

### **Nsfl1c (NM\_031981) Rat Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Nsfl1c (NM_031981) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Nsfl1c
Synonyms:	p47; XY40
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_031981
ORF Size:	1110 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200995).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_031981.2</a> , <a href="#">NP_114187.1</a>
RefSeq Size:	1450 bp
RefSeq ORF:	1113 bp
Locus ID:	83809
UniProt ID:	<a href="#">O35987</a>
Cytogenetics:	3q41
Gene Summary:	forms a complex with p97 that facilitates phosphatidylethanolamine-mediated membrane fusion; involved in reorganization of Golgi cisternae after mitosis [RGD, Feb 2006]



[View online »](#)