

Product datasheet for **RR200999L4V**

Cx3cl1 (NM_134455) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cx3cl1 (NM_134455) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Cx3cl1
Synonyms:	Cx3c; Scyd1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_134455
ORF Size:	1179 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200999).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_134455.1 , NP_604450.1
RefSeq Size:	3044 bp
RefSeq ORF:	1182 bp
Locus ID:	89808
UniProt ID:	O55145
Cytogenetics:	19p13
Gene Summary:	chemokine that binds its receptor CX3CR1; induces chemotaxis and increased intracellular calcium levels in microglia [RGD, Feb 2006]



[View online »](#)