

## Product datasheet for **RR201013L3V**

### **B4galt3 (NM\_001009539) Rat Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	B4galt3 (NM_001009539) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	B4galt3
Synonyms:	MGC72586
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001009539
ORF Size:	1185 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR201013).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001009539.1</a> , <a href="#">NP_001009539.1</a>
RefSeq Size:	2007 bp
RefSeq ORF:	1188 bp
Locus ID:	494342
UniProt ID:	<a href="#">Q6P768</a>
Cytogenetics:	13q24
Gene Summary:	Responsible for the synthesis of complex-type N-linked oligosaccharides in many glycoproteins as well as the carbohydrate moieties of glycolipids.[UniProtKB/Swiss-Prot Function]



[View online »](#)