

## Product datasheet for **RR201038**

### Ywhag (NM\_019376) Rat Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Ywhag (NM\_019376) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Ywhag  
**Synonyms:** 14-3-3G  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR201038 representing NM\_019376  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGTGGACCGGAGCAACTAGTGCAGAAAGCCCGCTGGCCGAGCAGGCGGAGCGCTACGACGATATGG  
CCGCGGCCATGAAGAACGTGACAGAGCTGAATGAACCACTGTCCAACGAGGAACGGAACCTCCTGTCGGT  
GGCCTACAAGAAGTGGTTGGGGCTCGCCGCTCCTCCTGGAGGGTCATCAGCAGCATCGAGCAGAAGACG  
TCTGCGGATGGCAACGAGAAGAAGATAGAGATGGTCCGAGCCTACCGGGAAAAGATTGAGAAGGAGCTGG  
AGGCCGTGTGCCAGGACGTGCTGAGCCTGCTGGACAACCTACCTGATCAAGAACTGCAGCGAGACCCAGTA  
CGAGAGCAAGGTGTTCTACCTGAAGATGAAAGGGGACTATTACCGCTACCTGGCCGAAGTGGCCACTGGG  
GAGAAAAGGGCGACCGTGGTGGAGTCGTCCGAGAAGGCCCTACAGCGAAGCCACGAGATCAGCAAGGAGC  
ACATGCAGCCCACCCACCCATCCGGCTGGGCCTGGCGCTCAACTACTCCGTTTTCTACTACGAGATCCA  
GAACGCCCCAGAGCAAGCCTGCCACCTGGCCAAGACCGCCTTCGATGACGCCATCGCCGAGCTCGACACT  
CTGAACGAGGACTCCTACAAGGACTCCACTCTGATCATGCAGCTGCTCCGAGACAACCTAACGCTCTGGA  
CGAGCGACCAGCAAGACGACGACGGCGGGAAGCAACAAC

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


**Protein Sequence:** >RR201038 representing NM\_019376  
Red=Cloning site Green=Tags(s)

MVDREQLVQKARLAEQAERYDDMAAMKNVTELENEPLSNEERNLLSVAYKNVVGARRSSWRVISSIEQKT  
 SADGNEKKIEMVRAYREKIEKELEAVCQDVLSLLDNYLIKNCSETQYESKVFYLLKMGDYRYLAEVATG  
 EKRA TVVESSEKAYSEAHEISKEHMQPHTPIRLGLALNYSVFYIEIQNAPEQACHLAKTAFDDAIAELDT  
 LNEDSYKDSTLIMQLLRDNLTLWTSQQDDDDGEGNN


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_019376  
**ORF Size:** 741 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_019376.2</a> , <a href="#">NP_062249.1</a>
<b>RefSeq Size:</b>	1634 bp
<b>RefSeq ORF:</b>	744 bp
<b>Locus ID:</b>	56010
<b>UniProt ID:</b>	<a href="#">P61983</a>
<b>Cytogenetics:</b>	12q12
<b>MW:</b>	28.3 kDa
<b>Gene Summary:</b>	belongs to the 14-3-3 family of proteins; mediates signal transduction via activation of protein kinase C and calcium/calmodulin-dependent protein kinase II [RGD, Feb 2006]