

Product datasheet for **RR201056L4V**

Akr1b1 (NM_012498) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Akr1b1 (NM_012498) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Akr1b1
Synonyms:	Akr1b3; Akr1b4; Aldr1; ALDRED; Alr; ALR-P-I; RATALDRED
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_012498
ORF Size:	948 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR201056).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_012498.1 , NP_036630.1
RefSeq Size:	1339 bp
RefSeq ORF:	951 bp
Locus ID:	24192
UniProt ID:	P07943
Cytogenetics:	4q22
Gene Summary:	regulates NF-kappa B mediated mitogenic signaling; may play a role in the polyol pathway [RGD, Feb 2006]



[View online »](#)