

## Product datasheet for **RR201081L3V**

### Hsd11b1 (NM\_017080) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Hsd11b1 (NM_017080) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Hsd11b1
Synonyms:	LRRGT00065
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_017080
ORF Size:	864 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR201081).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_017080.2</a> , <a href="#">NP_058776.2</a>
RefSeq Size:	1585 bp
RefSeq ORF:	867 bp
Locus ID:	25116
UniProt ID:	<a href="#">P16232</a>
Cytogenetics:	13q27
Gene Summary:	catalyzes the interconversion of cortisol and cortisone; plays a role in glucocorticoid metabolism; may regulate blood pressure [RGD, Feb 2006]


[View online »](#)