

Product datasheet for **RR201087**

Tbcc (NM_001108200) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Tbcc (NM_001108200) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Tbcc
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR201087 representing NM_001108200
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGGGAGTGGACTGCTCGGTGGCTCTGGCGGACGCTGCCCGGGTCCCCCGGGACCTGAGCCTGG
TACCCGAGCGACTTCAGAGGCGCGAGCAGGAGCGACAGATGGAGGTGGAGAGCGCAAGCAGAAGCGGCA
GGACCAGGAGGTGGAGGAGGAGAAGAGCGGCTTCTTCGCCCGCCCTTCGCCCGGAGCGCGCGGAGCG
GAGGAGCTGCTGCAGGGCGAGGCGTCCGCCGAGCGGCTGGAAGAGGCGGCAACCGGCTCCAGGCGCTGC
GGAAACTCTTGAACGACTCGGTTCTGTTCTCGCTGCCTACGACCTGCGGACGGCCAGGCGGCGTGGC
GCAGCTGCAGTCGGTCTAACGGAACGGCGCCAGGAGCTGCAGCCCAAGAAGCGCTTCGCTTCAAGGCT
CGAAAGAAGGATGCTGCCGGGACTGCCAGGTAGATGCCGCACCCGTCACCTCGGCGGCCCTTCTCCTC
CCTTGACCGAGGAGGAGGAGCTCCCGGGCCAGCTGGGCCTGTGGCTTCTCCAACCTGGAATCCAAGA
TCTGGAGAAGAGAGCGGAGGAGCTGCACCAGCGTGACATCTTTTGGAGCGACCTGACCAACTGCACTGTC
AGACTGCGTGGCAATCCAACACCTTGGGTTGGCAAAGGCCCGAGGCTGCAAGGTGCTCTGCGGCCAG
TGACCACCTCTGTTCCTGGAGGACTGTGAGATTGCATCCTGGCTGTGGCTTGCAGCAGCTCCGCGT
ACACACTACCAAAGACACCCGCATCTTCTTGCAGGTGACCAGCAGGGCCATCGTGGAGGACTGCAGCGGC
ATCCAGTTTGCTCCTTACACCTGGAGCTACCCTGGGATCGACAAGGACTTCCAGGACTCCGGTTTAGATA
GGAGCAAAAATAACTGGAACCAGGTTGACGATTCAATTGGCTAGCCAAAATGTGCCCTCCCAAAATTG
GAGTATTCTTCTGAGGAGGATCGAGACATCCAGTGGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR201087 representing NM_001108200
 Red=Cloning site Green=Tags(s)

MEGVDCSVALADAAAGSPRDL SLVPERLQRREQERQMEVERRKQKRDQDEVEEEKSGFFAAAFARERAAA
 EELLQGEASAERLEEAANRLQALRKL LND SVLFLAAYDLRQGQAALQLQSVL TERRQELQPKRF AFKA
 RKKDAAGTAQVDAAPVTSAPSPPLTEEGAPGASWACGFSNLESQDLEKRAEELHQRDILLSDL TNCTV
 RLRGNPNTLR LAKARGCKVLCGPVTT SVFLEDCRDCILAVACQQLRVHTTKDTRIFLQVTSRAIVEDCSG
 IQFAPYTWSYPGIDKDFQDSGLDRSKNNWNQVDDFNWLAQNVSPNWSILPEEDRDIQWD


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001108200

ORF Size:	1020 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001108200.2 , NP_001101670.1
RefSeq Size:	1184 bp
RefSeq ORF:	1023 bp
Locus ID:	316221
Cytogenetics:	9q12
MW:	38.1 kDa