

Product datasheet for RR201106L3

Kng111 (Kng1) (NM_001009628) Rat Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Kng111 (Kng1) (NM_001009628) Rat Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Kng1
Synonyms:	kininogen; Kng111
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR201106).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.

ACCN:	NM_001009628
ORF Size:	1290 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001009628.1 , NP_001009628.1
RefSeq Size:	1519 bp
RefSeq ORF:	1293 bp
Locus ID:	288001
UniProt ID:	P08932
Cytogenetics:	11q23
Gene Summary:	Kininogens are plasma glycoproteins with a number of functions: (1) as precursor of the active peptide bradykinin they effect smooth muscle contraction, induction of hypotension and increase of vascular permeability. (2) They play a role in blood coagulation by helping to position optimally prekallikrein and factor XI next to factor XII. (3) They are inhibitor of thiol proteases.[UniProtKB/Swiss-Prot Function]

