

## Product datasheet for RR201151L4

### Cask (NM\_022184) Rat Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Cask (NM_022184) Rat Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Cask
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR201151).
Restriction Sites:	Sgfl-RsrII
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.



[View online »](#)



RefSeq ORF: 2730 bp  
Locus ID: 29647  
UniProt ID: [Q62915](#)  
Cytogenetics: Xq12  
Gene Summary: may play a role in intracellular Ca(2+) mediated signaling pathways [RGD, Feb 2006]