

Product datasheet for RR201181L4

Ppif (NM_172243) Rat Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ppif (NM_172243) Rat Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Ppif
Synonyms:	CyP-D; CypD; PPIase
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR201181).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.



[View online »](#)

Locus ID: 282819
UniProt ID: [P29117](#)
Cytogenetics: 16p16
Gene Summary: mitochondrial isozyme of peptidyl-prolyl cis-trans-isomerase; enzyme activity is inhibited by cyclosporin A [RGD, Feb 2006]