

## Product datasheet for **RR201224L3V**

### Wif1 (NM\_053738) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Wif1 (NM_053738) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Wif1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_053738
ORF Size:	1095 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR201224).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_053738.1</a> , <a href="#">NP_446190.1</a>
RefSeq Size:	2182 bp
RefSeq ORF:	1098 bp
Locus ID:	114557
UniProt ID:	<a href="#">Q6IN38</a>
Cytogenetics:	7q22
Gene Summary:	human homolog has decreased mRNA expression in some prostate cancers and other cancers; may act as a secreted inhibitor of the Wnt signaling pathway [RGD, Feb 2006]



[View online »](#)