

Product datasheet for RR201253

Chka (NM_017127) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Chka (NM_017127) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Chka
Synonyms:	Chk; CK-R
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR201253 representing NM_017127 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAAACCAAGTTCTGCACCGGGGCGAGGCCGAGCCGTCCCCGCTTGGGCTGCTGCTGAGCTGCGGTG
GCAGCGCTGCCCGACGCCCGCGTAGGGCAGCAGCGCATGCCCGAGGGCAGCTGGAGTCCAAGCAGT
TGGTGGCCGGTCCCAACCTCTCGCGCTGCCGCCGCCACCACCGCCGCCCTGCCGCTGCCCGCCGCCA
TCACCGCCGCTAGCGGACGAACAACCCGAGCCCGGACGCGGCGCAGGGCCTACCTGTGGTGAAGGAAT
TCCTGCCCGGAGCCTGGAGGGGCCTTCGCGAGGACCAGTTCACATCAGTGTATCAGGGGTGGTCTCAG
TAACATGCTGTTCCAGTGTCTTCCAGACTCCATAGCCAGTGTGGTGTGAACTCGGAAAGTGCTC
TTGCGACTGTATGGGGCAATCTTAAAGATGGGGCTGAAGCAATGGTCTGGAGAGTGTATGTTTGCCA
TTCTTGACAGAGGTCACCTGGGCCAAAACCTCTATGGCATCTTCCGCAAGGCCGACTGGAGCAGTTTAT
CCCGAGCCGGCGATTGGACACTGAAGAATTATGTTTACCAGATATTTCTGCAGAAATAGCTGAAAAAATG
GCCACATTTATGGTATGAAAATGCCATTCAATAAGGAACCAAAATGGCTTTTGGAAACAATGAAAAAT
ACCTGAATCAAGTACTAAGACTTAAATTCAGCAGGGAGGCCAGAGTTCAACAAGTGCACAAGTTCCTCTC
TTACAATCTGCCTCTCGAGCTTGAGAACCTGAGGTCATTGCTGCAGTATACTAGATCCCCAGTTGTGTTT
TGTCATAATGACTGTCAAGAAGTAATATCTTATTGTTGGAAGGCCAAGAGAATTCTGAAAAGCAGAAGT
TGATGCTCATTGACTTTGAATACAGCAGTTACAATTACAGGGGATTTGACATTGGAATCATTTCTGTGA
ATGGATGTATGATTATACCTATGAAAAGTATCCTTTCTCAGAGCAAACATTCAGAAGTATCCTACCCGA
AAACAACAGCTCCATTTTATTTCAAGTTACTTGACTACATTCCAAAATGATTTTGAAGCCTCAGCAGTG
AAGAGCAGTCTGTACAAAAGAAGACATGTTGCTTGAAGTCAACAGATTTGCCCTGCCTCTCATTTCTC
CTGGGGACTTTGGTCCATTGTACAGGCCAAGATCTCATCCATTGAATTTGGGTACATGGAATATGCCAA
GCCAGTTTCGATGCTACTTTGACCAGAAGAGGAAGCTTGGTGTG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR201253 representing NM_017127
 Red=Cloning site Green=Tags(s)


MKTKFCTGGEAEPSPGLLLLSCGGSAAPTPGVGQQRDAAGELESKQLGGRSQPLALPPPPPPPLPLPPPP
 SPPLADEQPEPRTRRRAYLWCKEFLPGAWRGLREDQFHISVIRGGLSNMLFQCCLPDSIASVGDPRKVL
 LRLYGAILKMGAEAMVLESVMFAILAERSLGPPLYGIFPQGRLEQFIPSRRLDTEELCLPDISAEIAEKM
 ATFHGMKMPFNKEPKWLFGTMEKYLNQVLRKFSREARVQQLHKFLSYNLPLELENLRSLQLQYTRSPVVF
 CHNDCQEGNILLLEGQENSEKQKLMIDFEYSSYNYRGFDIGNHFCWEMYDYTYEKYPFFRANIQKYPTR
 KQQLHFISYLLTTFQNFESLSSEEQSATKEDMLLEVNRFFALASHFLWGLWSIVQAKISSIEFGYMEYQA
 ARFDAYFDQKRKLGV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN:

NM_017127

ORF Size:	1305 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_017127.1 , NP_058823.1
RefSeq Size:	2540 bp
RefSeq ORF:	1308 bp
Locus ID:	29194
UniProt ID:	Q01134
Cytogenetics:	1q43
MW:	49.7 kDa
Gene Summary:	kinase of the CDP-choline pathway [RGD, Feb 2006]