

Product datasheet for **RR201276L3V**

II17a (NM_001106897) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	II17a (NM_001106897) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	II17a
Synonyms:	CTLA-8; IL-17; IL-17A; II17
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001106897
ORF Size:	474 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR201276).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001106897.1 , NP_001100367.1
RefSeq Size:	1144 bp
RefSeq ORF:	477 bp
Locus ID:	301289
Cytogenetics:	9q13
Gene Summary:	cytokine; may be involved in immune system function or to cell death and cell survival [RGD, Feb 2006]



[View online »](#)