

Product datasheet for **RR201292L4V**

Rpap3 (NM_001004243) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rpap3 (NM_001004243) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Rpap3
Synonyms:	MGC94954
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001004243
ORF Size:	1977 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR201292).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001004243.1 , NP_001004243.1
RefSeq Size:	2220 bp
RefSeq ORF:	1980 bp
Locus ID:	300189
UniProt ID:	Q68FQ7
Cytogenetics:	7q36
Gene Summary:	Forms an interface between the RNA polymerase II enzyme and chaperone/scaffolding protein, suggesting that it is required to connect RNA polymerase II to regulators of protein complex formation.[UniProtKB/Swiss-Prot Function]



[View online »](#)